

NEW

NALGENE®

1L Centrifuge Bottle with Sealing Cap



Process a full liter at super speed

NALGENE 1L Centrifuge Bottle with Sealing Cap Assembly, polycarbonate or polypropylene bottle, polypropylene plug, modified PPE closure, and silicone gasket

- Leakproof up to 15,810 x G
- Excellent clarity and mechanical strength
- USP Class VI, non-cytotoxic
- Designed to fit F6S 6x1000Y, F8S 6x1000Y, and F9S 4x1000Y FIBER*Lite*® carbon fiber rotors

1L Centrifuge Bottle with Sealing Cap



Designed to boost productivity and pellet recovery in high speed, large volume gross separation applications. Process one full liter in each bottle at speeds up to 15,810 x G. Graduated, wide-mouth bottle simplifies decanting, re-suspensions and pellet removal. Smooth inner walls are less likely to shear delicate cells. Ideal for bacterial, yeast and tissue isolations. Improve viral harvesting in vaccine production. Reduce shearing in mammalian culture separations. Ergonomic design facilitates loading, recovery and storage of samples. Bottles are a perfect fit to FIBERLite F6S 6x1000Y, F8S 6x1000Y, and F9S 4x1000Y carbon fiber rotors. Recommended operating temperature range is from 4°C to 22°C; storage to -70°C. Autoclavable.



Super Speed Centrifuge Bottles

| NALGENE Cat. No. | Description | Nom. Cap., mL | Qty/Case |
|------------------|---|---------------|----------|
| 3140-1006 | 1L Super Speed Centrifuge Bottle, SCA, PC* | 1,000 | 6 |
| 3141-1006 | 1L Super Speed Centrifuge Bottle, SCA, PPCO** | 1,000 | 6 |

*NOTE: Use only Cat. Nos. DS3132-BLUE and 713140-00PC replacement sealing closure assemblies and o-rings with the 1L polycarbonate bottles Cat. No. 3140-1006.

**NOTE: Use only Cat. Nos. DS3132-ORED and 713141-PPCO replacement sealing closure assemblies and o-rings with the 1L polypropylene bottles Cat. No. 3141-1006.

Replacement SCA and o-rings

| NALGENE Cat. No. | Description | Size | Qty/Case |
|------------------|--|------|----------|
| DS3132-BLUE | Sealing Closure Assembly for 3140-100x ONLY | 85mm | 6 |
| 713140-00PC | Replacement BLUE O-rings, for 3140-100x ONLY | 85mm | 6 |
| DS3132-ORED | Sealing Closure Assembly for 3141-100x ONLY | 85mm | 6 |
| 713141-PPCO | Replacement RED O-rings, for 3141-100x ONLY | 85mm | 6 |



NALGENE®

www.nalgene.com

North America • Tel: 1-800-625-4327 • technical.nalgene@thermofisher.com

Asia Pacific • Tel: +65 6770 2807 • intlmtg@thermofisher.com

China • Tel: 86-21-68654588 • info.nnichina@thermofisher.com

Europe (NALGENE) • Tel: +44 (0) 1432 263933 • sales@nalgene.co.uk

Europe (NUNC) • Tel: +45 4631 2000 • info.nunc@thermofisher.com

India • Tel: +91-22-67162200 • contact.LPG.in@thermofisher.com

Japan • Tel: +81 3 3816 3355 • info@nalgenunc.co.jp • www.nalgenunc.co.jp

All other locations • (USA, International Department) Tel: +1 585 899 7198 • intlmtg@thermofisher.com



ISO 13485:2003
FM 31464